

B&G introduces a storage canister for the AccuSpray Elite. No need to refill your AccuSpray Elite while in the field. Just carry an extra AccuSpray canister that's already premixed with your chemical. Change-out is simple! Just remove the empty canister and replace with a full one. You'll be back performing treatments in less than a minute.

AccuSpray Elite Canister

- Stainless steel bottle
- Precision machined PVC cap
- Positive sealing cap insures leak free performance
- Heavy duty construction
- Great to use with the B&G IPM Bag



AccuSpray Holster

- Chemical and wear resistant material
- Belt loop for easy access
- Pockets to conveniently hold straws
- Heavy duty construction and economically priced



NUMBER	DESCRIPTION
24000115	AccuSpray Elite Spare Canister and Cap
24000072	Holster for AccuSpray Elite
45000132	IPM Bag with Optional Kit

For more information contact salesdept@bgequip.com or visit www.bgequip.com